ERKA.

Bedienungsanleitung
Instruction manual
Mode d'emploi
Instrucciones de manejo
Istruzioni per l'uso

Aneroid – Blutdruckmeßgeräte
Aneroid – Sphygmomanometers
Aneroid – Tensiomètres
Aneroid – Esfigmomanómetros
Aneroid – Sfigmomanometri

CE 0124

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Dear Customer,

With this blood pressure gauge you have chosen a quality product from ERKA.

We are delighted that you have shown your confidence in ERKA. quality and we thank you.

ERKA. instruments meet the valid European standard DIN EN 1060 and are in accordance with Annex V, Section 3, of the Directive 93/42 EEC dated June 14, 1993. This is certified by the CE-mark on the scale of the instrument. Each blood pressure gauge has a serial number and is registered accordingly.

The instrument is guaranteed against any material defects which may occur. The guarantee doesn't cover improper use, natural wear and tear or outside interference.

Please ensure that the card is stamped by your dealer.

2. Description of instrument

ERKA. aneroids work according to the auscultatory measuring method Riva Rocci / Korotkoff.

The necessary pressure in the cuff is created by the inflation bulb. The cuff is inflated by a few pumpings. The air release valve is sensitively adjustable and ensures a controlled pressure reduction during the measuring.

The reading scale has a measuring range from 0 - 300 mmHg with a graduation of 2 mmHg. The measuring accuracy over the total measuring range is ± 3 mmHg.

After calibration and a 100 % final control each manometer is sealed against unauthorised opening.

3. Preparations for measuring

Important ! Use the correct cuff-size !

As standard the instrument is supplied with an adult cuff size 13 x 55 cm. This cuff is suitable for upper arms with a circumference of 29 - 40 cm. For a larger or smaller circumference the corresponding special sized cuffs should be used. (Please see chapter 5.)

Fixing the cuff (general hint)

The cuff must be completely vented before being fixed. It should be tightly fixed to the left upper arm approx. 2 cm above the elbow joint. The center mark (black triangle on the cuff) must lie on the position of the arteria brachialis.

The arm should be positioned so that the upper part of the cuff is level with the aorta valve.

Fixing the cuff of home-care instruments (special hint)

The inflator tubes must point towards the hand. The microphone part of the stethoscope which is integrated in the cuff must be placed on the arteria. The free end of the cuff must be drawn through the metal bow and afterwards be closed by pressing it onto the velcro closure. After that the stethoscope has to be connected to the microphone part. With this model it is the microphone part and not the center mark which is placed on the arteria brachialis !

a) manometers with screwable air release valve

Take the instrument into your right hand and check that the air release valve is closed. (Turn clockwise between thumb and index finger.)

The cuff is inflated by pumping the bulb. The pointer of the instrument goes up and shows the respective pressure in the cuff. Increase the pressure to approx. 20 -30 mmHg above the expected systolic value. If you don't know this value, inflate to 240 mm Hg or higher. Wait for 4-5 pulse beats, then open the air release valve a little. The recommended descending rate is 2-3 mmHg per second.

Observe the pointer going down. As soon as pulsation is audible in the stethoscope, read the upper (systolic) value on the manometer. On continued air release a number of pulsations of decreasing intensity will follow. The last audible pulsation gives the lower (diastolic) value.

b) manometer with automatic valve

The handling of this instrument is easier and more comfortable. By pumping the bulb the hexagonal button in the centre of regulating screw comes up. Measuring should be made as described in chapter 4a.

The air release speed is fixed once. For that turn the regulating screw with the arrow mark.

a little to the left = air escapes more slowly
a little to the right = air escapes faster

The regulating screw does only one complete turn. Don't overturn it beyond the right or left stop.

After having taken the blood pressure, press the button in the centre of the regulation screw for a quick and complete venting of the cuff. Press this button too to stop the measuring immediately.

5. Choice of cuff-size

The following list is a recommendation for the use of cuffs with different sizes of rubber bags.

Patient	Circumf. of arm
Baby	8 - 13 cm
Infant	13 - 20 cm
Child	19 - 29 cm
Adult	29 - 40 cm
Adult large	38 - 50 cm
Extra-large	50 - 64 cm

We offer a wide range of cuffs of various sizes, colours, materials and fastening systems.

Standard instruments and cuffs contain latex components

For people with allergic reactions to latex material we offer latexfree alternatives.

- * The cuff pressure must not exceed 300 mmHg
- * The measuring time should not exceed 2 minutes
- * A repeated measure is to be taken only after 2 min. rest
- * Protect the instrument from breaking and heavy vibration
- * Protect the instrument from dirt and moisture
- * Avoid damaging the rubber parts from sharp or pointed objects
- * Store the instrument at a temperature between 10 and 40°C and a relative air humidity below 85 %
- * Do not expose the instrument, in particular the rubber parts, to direct ultraviolet radiation
- * The cuff is to be used with a hygienic protective strip or it must be cleaned before it is used for the next patient.

7. Maintenance and care

a) Manometer

In principle the manometer is maintenance-free. Dust particles in the air can block the filter in the main body after a certain time. The bulb can be removed for cleaning or replacement.

b) Cuff

Before cleaning the cuff, the rubber bag must be removed. Covers with velcro fastening are to be folded before washing so that the parts of fleece and of hooks are well attached.

Hints for care

	Non-shrink popelin colour: green	100% cotton black
Washing Ironing position Dry cleaning Dryer	60° mild washing weak as delicate clothing mild	60° standard normal as norm. clothing not possible
	easy to wipe PUR-laminated material all colours	
Wet wiping Dry cleaning Disinfecting	with a damp cloth,if necessary using a mild washing liquid as delicate clothing The cuffs have been tested to withstand the following recommended disinfectants: Cidex, Sporiciden, Microzid, Isopropyl-alcohol 70%, Ethanol 70%, Buratan liquid. After washing in disinfectants the cuff cover should be rinsed and allowed to air dry.	
Sterilization	The cuff has been tested to withstand Ethylene Oxide (EtO) gas sterilisation. Be sure that all safety precautions regarding aeration after EtO exposure are followed.	
	!!! Never autoclave !!!	
	Nylon / all colours	
Washing Wiping	30° mild washing + mild detergent with sponge or soft brush + mild detergent, rinse with water	

Blood pressure gauges must be calibrated and marked CE when they are put on the market.

We recommend for the doctor's and the patient's security that the instrument be checked and recalibrated at least every 2 years either by the manufacturer or by an authorised workshop.

The date of manufacture is shown by the first two numbers of the serial number.

Instruments which deviate from zero position in the range of 0 - 4 mmHg should be subject to immediate recalibration. Instruments without seal (for example after a repair) are considered as uncalibrated and have to be recalibrated.

We emphasize that only original ERKA. parts should be used for any spare or repair. A combination or attachment of ERKA. products with other makes is not possible.

We trust that working with your ERKA. blood pressure gauge will always meet with your satisfaction.

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