

P.M.S (Instruments) Ltd Tel 01628 773233 Fax 01628 770562
Waldeck House, Waldeck Road, Maidenhead, SL6 8BR

[Www.pmsinstruments.co.uk](http://www.pmsinstruments.co.uk)

Price list 2009 **cardiotest line**



ECG Machines

Software and Systems



Contents

Basic, Interpretive and Power ECGs

PC Compatible ECGs and Software

PC Link ECG

Accessories & Consumables

Excellence in a heartbeat – ensure quality traces from start to finish.

Based on experience and close co-operation with healthcare professionals and cardiology specialists we have ensured that the seca range of ECG machines comprehensively covers cardiology requirements in both primary and secondary care. The policy of continuous product improvement sees the introduction of model and software enhancements. The attractively priced CT8000i offers adult and paediatric interpretive software as a standard feature as well as computer connectivity and A4 print options. Software developments enable PC connection for archiving ECG records for paperless data management. seca Passport (Archimed) software will take the ECG and attach it within the individual patient's record enabling not just an on-screen image to be captured and viewed but also a choice of formats and measurements, etc. The file can be sent by e-mail as an image or sent as a file and by utilising an ECG viewer, ensures confidentiality. With increased connectivity comes user choice; as an example printers can be used to provide inexpensive A4 formats where preferred.

CT3000i Interpretive

- 12 lead, 3 Channel ECG with high quality interpretive programme.
- Auto-start button automatically produces the ECG traces complete with Rhythm Strip, Measurements and Interpretation.
- Filters for baseline, electrical and muscle tremor interference, providing clear traces.
- Lead off warning and identification.
- Useful copy button for extra ECG traces.
- 12 simultaneously acquired ECG traces printed on 90 mm paper, eliminates trace overlap.
- Paper tray accommodates z-fold or roll paper.
- Manual button provides real time print-outs.
- Mains and battery operation.
- Battery recharges automatically and when fully charged provides 4 hours of continuous operation.




By pushing the auto-start button the ECG traces are produced automatically.

- Secure and portable within the seca hard carry case option (see accessories and consumables).
- Supplied complete with mains cable, patient cable, pack of z-fold paper, 500 disposable electrodes and alligator clips.
- Dimensions: 293 x 210 x 73 mm
- Weight: 2.9 kg

Price

CT3000i interpretive £ 1,275

High quality print out.



CT6i Interpretive

A4 Presentation.

- 12 Channel ECG with high quality interpretive programme.
- Choice of up to 12 leads per page or the popular 4 x 3 leads plus Rhythm Strip, Intervals and Interpretation on one A4 page.
- Auto-start one touch operation or manual real time.
- Large A4 paper tray providing up to 240 ECG pages.
- Mains and battery operation with automatic recharge.
- Fully charged battery will provide up to 400 ECGs.
- 12 lead simultaneous acquisition with lead off warning and identification.
- Interference filters for baseline drift, electrical interference and muscle tremor.
- Copy button for extra traces.
- Supplied complete with mains cable, patient cable, pack of z-fold paper, 500 disposable electrodes and alligator clips.
- Dimensions: 399 x 338 x 104 mm
- Weight: 5 kg



seca CT500 ECG Trolley

Price

CT6i interpretive £ 1,975




CT6Pi Interpretive

A4 Presentation.

- A4 size 12 channel ECG with high quality adult and paediatric interpretive programme. Inbuilt RS232 port for computer connection.
- LCD screen for viewing traces, stored ECGs and input of patient details via QWERTY keyboard.
- Auto-start button automatically produces the ECG traces complete with rhythm strip, measurements and interpretation on one A4 page (or other formats as required).
- 12 lead simultaneous acquisition with interference filters for base line drift, electrical and muscle tremor.
- Mains and battery operation. Battery recharges automatically, fully charged will provide up to 3 1/2 hours continuous operation.
- On-screen lead off warning and Identification.
- Copy button for extra traces.
- Supplied complete with mains cable, patient cable, pack of paper sufficient for 240 ECGs, 500 disposable electrodes and alligator clips.
- Dimensions: 399 x 338 x 104mm
- Weight: 5 kg


Price
CT6Pi £ 2,750



Optional: paperless ECG records.



Patient data input screen.



Traces and lead off indicator.

CT8000i Interpretive

A4 print Connection.

Compact, light weight and multifunctional ECG offering the best of both worlds –

- Paperless ECG option – download and electronically store ECG traces inside individual patients' records via RS232 port.
- Stand alone use and PC connectivity for paperless ECG options.
- 12 lead, 3 Channel high quality interpretive ECG with LCD screen and alpha-numeric keyboard.
- Choice of both internal and external print formats, including automatic rhythm strip.
- One touch automatic or instant real time operation.
- Internal memory for up to 45 stored ECG records.
- High quality adult and paediatric interpretive software.
- Digital filters for mains, base line and muscle tremor.
- Easy-to-use one-touch button operation.
- Additional advanced filters (SSF and SBS) for clear traces.




- Dedicated soft keys for ease of storing, printing and transmitting.
- Integral port to connect to a standard DeskJet or Laser printer for an A4 print out.
- Optional communication module (built-in Lan and modem connections).
- Mains/rechargeable battery operation.
- Secure and portable within the seca hard carry case option (see accessories and consumables).
- Supplied complete with mains cable, patient cable, pack of z -fold paper, 500 disposable electrodes and alligator clips.
- Dimensions: 285 x 195 x 70 mm
- Weight: 2.5 kg

Price
CT8000i Interpretive £ 1,795
Thrombolysis Upgrade £ 265



High quality print out.



Optional: paperless ECG records.

CT8000P Interpretive

A4 Presentation.

- Equipped to provide paperless ECG options when combined with seca Archimed Passport Software.
- 12 Channel high quality interpretive ECG with LCD screen and alpha-numeric keyboard.
- High quality adult and paediatric interpretive software.
- Easy selection of traces from 12 x 1, 6 x 2 or the popular 4 x 3 plus Rhythm Strip one page report.
- On-screen lead test and lead-off indication.
- One touch automatic or instant real time operation.
- High quality digital filters for clear traces (mains, base line and muscle tremor).
- Internal memory for up to 45 stored ECG records.
- All parameters including filters are easily programmed through the on screen menu for a one touch ECG.
- SVGA output for external colour monitor connection.



Dedicated one-touch keys.



Powered papertray.



- Dedicated one touch keys for storing, printing and transmitting.
- Easy-load paper tray (220 A4 pages).
- Upgradeable Thrombolysis Indication, please enquire for full details.
- Optional communication module (built-in Lan and modem connections).
- Mains/rechargeable battery operation.
- Supplied complete with mains cable, patient cable, pack of z-fold paper, 500 disposable electrodes and alligator clips.
- Dimensions: 330 x 381 x 76 mm
- Weight: 5 kg

Price

CT8000P Interpretive	£ 2,995
Thrombolysis Upgrade	£ 265



CT8000L Interpretive

A4 Presentation.

Offering all the features of the CT8000P excluding on-screen tracing preview.

- Equipped to provide paperless ECG options when combined with seca Archimed Passport Patient Record Software.
- 12 Channel high quality interpretive ECG with LCD screen for menu navigation and alpha-numeric keyboard.
- High quality adult and paediatric interpretive software.
- Easy selection of traces from 12 x 1, 6 x 2 or the popular 4 x 3 plus Rhythm Strip one page report.
- On-screen lead test and lead-off indication.
- One touch automatic or instant real time operation.
- High quality digital filters for clear traces (mains, base line and muscle tremor).
- Internal memory for up to 45 stored ECG records.
- All parameters including filters are easily programmed through the on-screen menu for a one touch ECG.
- Dedicated one touch keys for storing, printing and transmitting.

- Optional Communication module (built-in Lan and Modem connections).
- Easy-Load paper tray (220 A4 pages).
- Upgradeable option for Thrombolysis indication and/or full tracings preview.
- Mains/rechargeable battery operation.
- Supplied complete with mains cable, patient cable, pack of z-fold paper, 500 disposable electrodes and alligator clips.

Price

CT8000L Interpretive	£ 2,350
Thrombolysis Upgrade	£ 265





CT8000t Tele

New compact, light weight and multifunctional ECG offering unique Tele-medicine features:

- Stand alone use and PC connectivity for paperless ECG options.
- 12 lead, 3 Channel high quality interpretive ECG with LCD screen and alpha-numeric keyboard.
- Choice of both internal and external print formats, including automatic rhythm strip.
- One touch automatic or instant real time operation
- Internal memory for up to 45 stored ECG records.
- Tele-medicine applications available from the built in GSM module, which includes dedicated antenna and a slot for standard SIM card.
- Paperless ECG option – download and electronically store ECG traces inside individual patients' records via RS232 or LAN ports.
- High quality adult and paediatric interpretive software.
- Digital filters for mains, base line and muscle tremor.
- Easy-to-use one-touch button operation.
- Additional advanced filters (SSF and SBS) for clear traces.



- Dedicated soft keys for ease of storing, printing and transmitting.
- Integral port to connect to a standard Deskjet or Laser printer for an A4 print out.
- Mains/rechargeable battery operation.
- Secure and portable within the seca hard carry case option (see accessories and consumables).
- Supplied complete with mains cable, patient cable, pack of z-fold paper, 500 disposable electrodes and alligator clips.
- Package includes, organiser, software in a single server configuration. Network and simcard charges not included.

Price

CT8000t Tele £ 4,850

Software

Paperless ECG Records.

The seca CT6Pi and the CT8000 range of ECG machines are designed to print (or not if required), store to memory and transmit into patient record systems. The seca Archimed Passport Software enables not only the ECG trace, but also a choice of format and measurements etc to be attached into the individual patient record. Once in the computer system the trace can be printed, faxed or emailed as desired.

- Transmission of the ECG trace and patient data through one touch function key on the ECG machine.
- User selection of viewing choices (e.g. One page report or 3 x 4, etc.).
- All lead measurement figures included in patient record.

Price

Archimed Single User	£ 99
Archimed Passport Patient Record Package	£ 299



Archimed Single User

This software operates on one computer, archiving ECG data, and is suitable for the single handed practice.

Archimed Passport Patient Record Package

This software receives the trace from the ECG machine and enables attachment to the individual patient record for storage and viewing throughout the practice network.



CT6Pi Interpretive



CT8000i Interpretive



CT8000P Interpretive



CT8000L Interpretive



CT110 CardioConcept

Interpretive PC based ECG.

Using a conventional 10 lead ECG cable, the seca CT110 CardioConcept acquires the ECG from the patient directly into the treatment room computer.

- The interface is easy to use and the ECG recording can be viewed, printed and archived into the practice Patient Record system.
- High quality recordings (2000Hz sampling) of all twelve leads for viewing, printing or comparison with previous tracings.
- Complete with ECG interpretive software.
- Robust defibrillator proof connection to PC, the unit is powered from the USB port, eliminating the usual mains cable.
- Requiring only two clicks of the mouse to take an ECG recording.



Easyload software

- The seca CardioConcept interface is user friendly and simple to operate.
- Paperless ECG records. If hard copy is required, the ECG can be printed complete with grid onto plain copier paper.
- Archive software identifies and filters data for retrieval of specific ECG classifications.
- Data can be attached to all leading Patient Record systems. Data export choices inc. Jpeg, tif, PDF, etc..
- seca CT110 CardioConcept supplied complete with USB connecting cable, alligator clips, 500 disposable electrodes, full length traditional patient cable and CardioConcept software on USB flash drive (requires Windows 2000 or above).

Price
CT110 CardioConcept £ 1,995



seca CT500 ECG Trolley

- Dedicated for use with the seca range of ECG machines.
- Spring loaded arm for control of patient cable and storage basket for electrodes and paper.
- Non-metallic top.

Price
seca CT500 ECG Trolley £ 249



seca 582 Hard Carry Case.

Cardiotest Accessories and Consumables

Price		
seca 582		£ 69.00
hard aluminium carry case for CT3000 range and CT8000i		
seca 580		£ 39.00
carry case for CT3000 range and CT8000i		
seca 426		£ 30.00
carry case for CT6 range and CT8000P		
seca Q-trace	Pack of 500	£ 29.00
disposable electrodes		
seca Alligator Clips	per pack	£ 18.50
4mm for disposable electrodes (pack of 10)		
Snapclip Connectors		£ 29.00
for disposable snap and tab electrodes (fits 3mm and 4mm patient cables)		
Suction Cup Chest Electrodes	each	£ 8.50
reusable (replacements)		
Limb Clamp Electrodes	each	£ 8.50
reusable (replacements) - single		
Limb Clamp Electrodes	per pack	£ 32.00
reusable (replacements) - colour coded (pack of 4)		
ECG Cream	each	£ 2.00
for use with reusable electrodes		
ECG Gel	each	£ 2.00
for use with reusable electrodes		



seca Q-Trace and Snapclip Connectors.



seca Q-Trace

ECG Paper

Price		
CT446	each	£ 10.00
CT3300, CT3600 - z-fold		
CT448	each	£ 10.00
CT6300, CT6600 - z-fold		
CT467	Pack of 10	£ 16.00
CT1000 - roll		
CT463R	each	£ 3.20
CT3000 range - roll (width 90 mm)		
CT463Z	each	£ 3.60
CT3000 range - z-fold (width 90 mm)		
CT480Z	each	£ 3.60
CT8000i - z-fold		
CT486Z	each	£ 15.00
CT6 range - z-fold (A4)		
CT480ZPi	each	£ 15.00
CT8000P - z-fold (A4)		

Minimum order value on consumables and paper £ 15.00.

The right to make modifications is reserved. These modifications always bring improvement.



All seca ECG machines are CE marked and meet all the legal requirements and european directives for such instruments.

3 year warranty applies to all models in the CT3000, CT6 and CT8000 ranges. 2 year warranty for all other models.

Contact your usual supplier and specify which seca ECG machine you wish to purchase or have demonstrated.

Valid from 1st March 2009.

Former prices are hereby cancelled. Please state complete machine specification when ordering including options as indicated.

Title: The title in any goods or accessories supplied will not pass to the buyer until payment has been made in full by the buyer to the Company in accordance with the terms of any contract to which the conditions apply.

Terms: Prices are net, including packing, but not VAT which will be added at the appropriate rate if applicable. Carriage/Delivery charges as necessary.



seca group

- seca deutschland
- seca france
- seca united kingdom
- seca north america
- seca schweiz
- seca zhong guo
- seca nihon
- seca mexico
- seca austria