

ERKA.SWITCH

**A new concept for blood pressure measurement in hospitals
Low cost combined with simple operation for better service
and control**



original size

- colour labels for individual marking
- easily adjustable for left or right hand use
- quick exchange of manometer housing
- calibrated cuff
- quick local recalibration
- 4 years guarantee on handle

ERKA. SWITCH

CE 0124

The new concept of ERKA. blood pressure measuring as introduced with the SWITCH is ideal for hospitals. Thanks to its excellent ergonomic design it is very comfortable to hold and it can be easily and rapidly adapted from right to left hand use. Moreover it is small and light enough to keep carrying with you all day long.



The new ERKA. SWITCH is a great design. It looks modern and by its colourful label I can easily identify which is mine. I particularly like the fact that we do not have to send them away for service and control – it used to take quite some time before we received them back.



The new ERKA. SWITCH is excellent value for money as it means I only have to buy the housing when service is needed. The idea is particularly appealing as there is no longer a sending back and forth involved. That not only saves money but takes also environmental issues into consideration.



With the new ERKA. SWITCH the official regular control is easy and can be done here in the hospital. If the equipment does not live up to the required standards I simply exchange the housing. That saves time and money. And polycarbonate is a highly shock resistant material which makes the SWITCH quite robust. It's tested according to AAMI SP 9 drop test standards and has a guarantee of 4 years on the hand part. Really a good thing for hospitals!



ERKA. is a family owned company with a tradition dating back over 115 years. From the very beginning ERKA. has largely contributed to the development of this product.

ERKA. manometers are excellent quality exclusively produced in Germany.

ERKA.

KALLMEYER Medizintechnik
GmbH & Co. KG
Im Farchet 15
D-83646 Bad Tölz
Rev.Std 001/Aug. 2004

phone: +49 8041 8009-0
fax: +49 8041 8009-39
E-Mail: switch@erka.org
Internet: www.erka.org

UK Distributor
PMS (Instruments) Ltd
Waldeck House
Waldeck Road
Maidenhead
Berks, SL6 8BR
Tel 01628 773233
Fax 01628 770562
web www.pmsinstruments.co.uk