

Type	Stethoscope with flat chest piece from aluminium, anodized for adults		
	UMDNS Code	13-755	
	Medical Product	Class 4	
Content of delivery	Stethoscope complete in decorative cardboardbox		
Dimensions of Stethoscope	length (total)	700	(mm)
	length of molded Y-tube	610	(mm)
	diameter of binaurals	5	(mm)
	length of binaurals	170	(mm)
	max. opening of binaurals	290	(mm)
Dimensions of Chest piece	diameter of chest piece	47	(mm)
	height of chest piece	10	(mm)
	diameter of diaphragm	44	(mm)
Spring power of binaurals	0,6 mm spring at 150 mm angle width	320	(g)
	at 200 mm angle width	450	(g)
Material	chest piece	aluminium, anodized	
	binaurals	brass, polished and chromium-plated	
	ear tips	PUR hard, colour black	
	diaphragm	PET, 0.25 mm thick	
	molded Y-tube	soft PVC	
		choice of 6 colours, latex free, resistant against common detergents and disinfection	
Weight	chest piece	25	(g)
	complete without packing	110	(g)
Dimensions of box	(L x W x H)	280 x 95 x 35	(mm)
Freight data	gross weight (incl. packing) / piece	200	(g)
	volume / piece	930	(cm ³)
Individual marking	of chest piece on request		
Warranty	2 years		
Manufacturer	ERKA. Kallmeyer Medizintechnik GmbH & Co KG Im Farchet 15 • D 83646 Bad Toelz Tel. +49 (0)8041 8009-0 • Fax: +49 (0)8041 8009-39 • E-Mail: info@erka.org		
Country of origin	Germany		