

Type	Stethoscope with double chest piece, switchable, from brass, chromium-plated, for adults		
	UMDNS Code	13-755	
	Medical Product	Class 4	
Content of delivery	Stethoscope complete with 2 spare ear tips in decorative cardboardbox		
Dimensions of Stethoscope	length (total)	690 (mm)	
	length of molded Y-tube	510 (mm)	
	diameter of binaurals	5 (mm)	
	length of binaurals	180 (mm)	
	max. opening of binaurals	300 (mm)	
Dimensions of Chest piece	diameter of chest piece	47 (mm)	
		bell shaped side	30 (mm)
	height of chest piece	22 (mm)	
	diameter of diaphragm	45 (mm)	
Spring power of binaurals	0,6 mm spring	at 150 mm angle width	260 (g)
		at 200 mm angle width	300 (g)
	0,7 mm spring	at 150 mm angle width	450 (g)
		at 200 mm angle width	600 (g)
Material	chest piece / binaurals	brass, polished and chromium-plated	
	ear tips super soft	silicone rubber, soft, colour light-grey	
	diaphragm	PET, 0.25 mm thick	
	non-chill-ring	neopren	
	molded Y-tube	soft PVC	
		colour black, latex free, resistant against common detergents and disinfection	
Weight	chest piece	93 (g)	
	complete without packing	169 (g)	
Dimensions of box	(L x W x H)	395 x 186 x 30 (mm)	
Freight data	gross weight (incl. packing) / piece	325 (g)	
	volume / piece	2200 (cm ³)	
Individual marking	of chest piece on request		
Warranty	3 years		
Manufacturer	ERKA. Kallmeyer Medizintechnik GmbH & Co KG Im Farchet 15 • D 83646 Bad Toelz Tel. +49 (0)8041 8009-0 • Fax: +49 (0)8041 8009-39 • E-Mail: info@erka.org		
Country of origin	Germany		