

HF [bpm]

60

Patient: James, Basil
Date of birth 01.01.1978 (33 Years)
Patient:
Sex/Size/Weight: unknown / -- cm / -- kg
Pacemaker:

Recording from: 08.03.2011 13:59:30
Blood pressure: -- / -- mmHg

P.M.S. (instruments) Ltd.
Waldeck House, Waldeck Rd.
SL6 8BR Maidenhead, Berks
Recorder
Finder




HF [bpm]

60

Patient: James, Basil
 Date of birth 01.01.1978 (33 Years)
 Patient:
 Sex/Size/Weight: unknown / -- cm / -- kg
 Pacemaker:

Recording from: 08.03.2011 13:59:30
 Blood pressure: -- / --mmHg

P.M.S. (instruments) Ltd.
 Waldeck House, Waldeck Rd.
 SL6 8BR Maidenhead, Berks
 Recorder
 Finder




Finding

rhythm
 sinus-rhythm
 specific findings
 Q in : V1 V2 V3
 broad R in : V1
 small S in : V1 V2
 QT-Dispersion : Std 0 ms Max 2 ms from 8/8 leads
 QT dispersion : Std 0 ms Max 2 ms from 6/6 CH-leads
 possible right atrial enlargement


QRS-T evaluation
 possible anterior myocardial infarction, probably old
 Summary
 possible abnormal ECG

HF 60 /Min.
 RR 1000 ms
 P 108 ms
 PQ 170 ms
 QRS 92 ms
 QT 388 ms
 QTc(B) 388 ms
 QTc(F) 388 ms
 QTr(H) 99 %



P 49°
 QRS 49°
 T 48°

	Duration [ms]			max. amplitude [mV]							[mV/s]		
	Q	R	S	P _{pos}	P _{neg}	Q	R	S	T _{pos}	T _{neg}	J	ST	ST
I	22	50	20	0.16		-0.08	0.77	-0.08	0.32				
II	22	50	20	0.24		-0.12	1.16	-0.12	0.48				
III	16	50	18	0.08		-0.04	0.38	-0.04	0.16				
aVR		22	50		-0.20		0.10	-0.96		-0.40			
aVL	14	48	10	0.04		-0.02	0.20	-0.02	0.08				
aVF	22	50	20	0.16		-0.08	0.77	-0.08	0.32				0.10
V1	18	48	20	0.08		-0.04	0.42	-0.04	0.18				
V2	22	50	20	0.17		-0.09	0.82	-0.09	0.34				
V3	22	50	20	0.24		-0.12	1.16	-0.12	0.48				
V4	22	50	20	0.34		-0.18	1.66	-0.16	0.69				0.10
V5	22	50	20	0.24		-0.12	1.17	-0.12	0.49				
V6	22	50	20	0.17		-0.09	0.83	-0.08	0.34				



HF [bpm]

--- **60**

--- mmHg

--- mmHg

Patient:

James, Basil

Date of birth

01.01.1978 (33 Years)

ECG Comparison

ECG recording: 08.03.2011 13:59:30



P.M.S. (instruments) Ltd.

Waldeck House, Waldeck Rd.


SL6 8BR Maidenhead, Berks

Recorder


Finder



Representative cycle - II



ST-Overview



Comparison ECG	Actual ECG
ST level(II): mV	mV
ST slope(II): mV/s	mV/s
RR: ms	1000ms
P: ms	108ms
PQ: ms	170ms
QRS: ms	92ms
QT: ms	388ms
QTc(B): ms	388ms
QTc(H): %	99%
QTc(F): ms	388ms

