



Medtronic
PHYSIO-CONTROL

Wall Cabinet

For use with the LIFEPAK® 500 AED

Having trained lay rescuers equipped with the LIFEPAK 500 AED placed conveniently throughout a facility can help save precious minutes and improve survival rates for cardiac arrest victims. The AED Wall Cabinet is ideal for facilities such as airports, health clubs, golf courses, hotels, high rise buildings, gated communities and sprawling manufacturing plants.

This cabinet offers three trim styles: semi-recessed (3" return trim), recessed (flat trim) and surface mounted. These three options allow for any mounting application. Cabinet is manufactured of steel finished with a baked on white epoxy finish.

The cabinet is equipped with a standard alarm that can be armed or disarmed with key lock. It can also be fitted with magnetic contacts and lead wires to tie into most new or existing low voltage fire and/or security systems. Standard labeling, acrylic glazing and pull handle are also included. Each cabinet includes one AED Wall Placard and CPR/AED Wall Poster.



Semi-recessed





Surface Mount



Recessed Flat Trim

Exploded top and side views refer to dimensions table



Model	Trim Style	TUB I.D.			Frame O.D.		Wall Opening		
		A	B	C	E	F	W	H	D
3012604-000	Flat Trim	14"	14"	6"	17 3/8"	17 3/8"	14 1/2"	14 1/2"	6"
3012604-001	3"	14"	14"	6"	17 3/8"	17 3/8"	14 1/2"	14 1/2"	3 1/8"
3012604-002	Surf. Mtd.	17 3/8"	17 3/8"	6"	17 3/8"	17 3/8"	-	-	-

The Americans with Disabilities Act (A.D.A.) of 1990 mandates that cabinets or other accessories not extend into corridor areas more than 4". If this applies to an existing building, model 3012604-000 or 3012604-001 must be used.

Includes AED Wall Placard and CPR/AED Wall Poster.



Medtronic Physio-Control
11811 Willows Road NE
P. O. Box 97006
Redmond, WA 98073-9706 USA
Tel: 425.867.4000
Fax: 425.867.4121
Internet: www.physiocontrol.com
Internet: www.medtronic.com

Europe
Basingstoke, Great Britain
Tel: 44.1256.782.727
Fax: 44.1256.782.728

Canada
Richmond Hill, Ontario
Tel: 905.709.4330
Fax: 905.709.3336

United Kingdom and Ireland
Basingstoke, Great Britain
Tel: 44.1256.782.727
Fax: 44.1256.782.728

France
Lyon - St Priest, France
Tel: 33.4.72.79.26.26
Fax: 33.4.72.79.26.36

Germany, Switzerland
Fernwald, Germany
Tel: 49.6404.91.41.0
Fax: 49.6404.91.41.20

Austria
Vienna, Austria
Tel: 43.1.240.44.160
Fax: 43.1.240.44.600

Italy
Sesto Fiorentino, Italy
Tel: 39.055.302.701
Fax: 39.055.302.4029

Netherlands
Hoofddorp, The Netherlands
Tel: 31.20.6.533.640
Fax: 31.20.6.535.822

Spain
Madrid, Spain
Tel: 34.1.375.6050
Fax: 34.1.375.6055

Scandinavia
Järfälla, Sweden
Tel: 46.8.580.945.00
Fax: 46.8.580.945.05

Asia Pacific
Christchurch, New Zealand
Tel: 64.3.3794.429
Fax: 64.3.3792.374

Latin America
Miami, Florida
Tel: 786.242.6661
Fax: 786.242.4505

Middle East
Dubai, UAE
Tel: 971.4.282.6532
Fax: 971.4.282.7970

Hungary
Budapest, Hungary
Tel: 36.1.214.2228
Fax: 36.1.214.2230

Poland
Warsaw, Poland
Tel: 48.22.613.38.13
Fax: 48.22.613.38.11

Czech Republic
Prague, Czech Republic
Tel: 420.2.2017.2277
Fax: 420.2.2056.1617

People's Republic
of China
Beijing, China
Tel: 86.10.6708.1162
Fax: 86.10.6708.1163