



P.M.S (Instruments) Ltd
 Waldeck House, Waldeck Road
 Maidenhead, Berkshire, SL6 8BR
 Tel: 01628 773 233 Fax: 01628 770 562

Patient data			
name	dipper ABPM-04		
ID	2345-6789	date of birth	
Mother name			
Age		Sex	

Verifies NICE quality criteria met

Quick identification of average daytime

Next report
 The ABPM was worn by the patient for 24 hours. During this period the average daytime blood pressure was 169/118 and the pulse rate was 88 beats per minute. According to NICE guidelines the patient exhibits stage 2 hypertension. The examination was divided into day 06:00-22:00, night 22:00-06:00. 64 successful readings were taken during the daytime period and the NICE key quality requirement for 14 has been met. During the daytime there were 64 readings (100.0%) above the NICE guideline figure of 135/85. The highest reading was 202/134 recorded at 9:15 22/05/1997. During the night time period 24 successful readings were taken and the average blood pressure was 135/87.


Normotensive, stage 1 or stage 2 hypertension

Automatic identification of highest reading and time

The reference for the NICE guideline is CG127 Hypertension: NICE guideline August 2011.

www.nice.org.uk/guidance/cg127

Colour trend graph



ABPM basic data (22/05/1997 09:00 - 23/05/1997 09:00)				
period	Overall	Morning	Day	Night
time range	Overall	06:00 - 06:00	06:00 - 22:00	22:00 - 06:00
SBP/DBP maximum limit value	~127/77~	-/-	135/85	120/70
SBP/DBP weighted average	158/108	-/-	169/118	135/87
Pulse weighted average	84	-	88	78
readings	88	-	64	24
successful measurement ratio %	99	-	98	100
SBP/DBP D/N index %	20/26			
morning surge	45			

Blood pressure list						
Time	T	SBP mmHg	DBP mmHg	pulse bpm	remark	information
09:00	A	189	131	077		
09:15	A	202	134	074		
09:30	A	201	128	085		
09:45	A	167	130	081		
10:00	A	181	117	085		
10:15	A	180	124	086		
10:30	A	178	127	079		
10:45	A	161	119	083		
11:00	A	180	123	091		
11:15	A	175	119	073		
11:30	A	190	123	094		
11:45	A	167	105	075		
12:00	A	154	105	070		
12:15	A	159	104	072		
12:30	A	184	131	087		
12:45	A	163	119	085		
13:00	A	190	125	084		
13:08	E					pressing of the event button
13:15	A	174	115	075		
13:30	A	173	116	074		
13:45	A	155	089	075		
14:00	A	144	100	079		
14:15	A	162	097	084		
14:30	A	137	101	080		
14:45	A	140	109	099		
15:00	A	180	130	099		
15:15	A	174	124	103		
15:30	A	181	124	109		
15:45	A	193	124	105		
16:00	A	170	113	095		
16:15	A	173	115	111		
16:30	A	160	115	105		
16:45	A	178	116	099		
17:00	A	168	114	108		
17:15	A	162	119	104		
17:30	A	175	122	095		
17:45	A	171	113	095		
18:00	A	166	112	094		
18:15	A	177	122	097		
18:30	A	160	114	093		
18:45	A	152	115	089		
19:00	A	169	125	097		
19:15	A	168	118	094		
19:30	A	158	117	088		
19:45	A	147	115	085		
20:00	E					pressing of the event button
20:00	A					cancelled by user
20:15	A	152	118	083		
20:30	A	162	119	072		
20:45	A	164	118	090		
21:00	A	150	112	084		
21:15	A	162	118	081		
21:30	A	164	116	091		

Patient name: **dipper ABPM-04**
 Patient ID: 2345-6789
 Date of birth:

ABPM report
 Start of the examination
 22/05/1997

Blood pressure list						
Time	T	SBP mmHg	DBP mmHg	pulse bpm	remark	information
21:45	A	158	116	079		
22:00	A	150	113	088		
22:20	A	169	110	080		
22:40	A	135	095	074		
23:00	A	123	084	077		
23:20	A	137	078	086		
23:40	A	121	074	081		
00:00	A	129	073	081		
00:20	A	122	079	083		
00:40	A	119	076	076		
01:00	A	134	097	096		
01:02	E					pressing of the event button
01:20	A	129	076	080		
01:40	A	145	084	078		
02:00	A	140	088	083		
02:20	A	140	093	073		
02:40	A	138	093	076		
03:00	A	142	095	079		
03:20	A	121	078	070		
03:40	A	137	083	069		
04:00	A	131	082	069		
04:20	A	127	088	080		
04:40	A	125	087	070		
05:00	A	137	085	073		
05:20	A	143	095	076		
05:40	A	165	108	075		
06:00	A	151	116	073		
06:15	A	172	114	077		
06:30	A	166	125	090		
06:45	A	161	127	085		
07:00	A	167	122	084		
07:15	A	181	134	087		
07:30	A	176	132	082		
07:45	A	190	119	086		
07:59	E					pressing of the event button
08:00	A	171	129	086		
08:15	A	170	124	084		
08:30	A	167	124	093		
08:45	A	181	123	087		
09:00	A	182	132	098		
09:52	E					pressing of the event button