

insight

Product Update

Winter 2011

Medset Flashlight PC Based ECG System

Medset Flashlight



Price
£1595.00

Was
~~£1795.00~~

- Powered by USB port. No batteries. Wireless Bluetooth version also available
- Uses validated German Hanover HES ECG interpretation protocol developed by Cardiologists to provide accurate ECG measurements and interpretation
- Generates high quality PDF or JPEG reports which can be attached to EMIS, SystemOne, INPS clinical management systems for a paperless solution
- Anti drift baseline, muscle and mains power filters for optimum high quality ECGs
- Comparison ECG feature
- Software checks ECG quality before recording
- Record up to 180 seconds of ECG
- 12 channel ECG system with user configurable reports

SB-220 Professional Spot Check Finger Oximeter

SB-220



Price
£77.00

Was
~~£99.00~~

- Bright OLED two colour display is easy to read
- Displays SPO2, heart rate and pulse strength
- Changeable display viewing direction
- 2 AAA batteries - last for up to 16 hours continuous use
- Audible alarm function if SPO2 below 90%
- Weighs just 37g size 63.5 x 34 X 35mm

Prices exclude carriage and VAT. Offers available until 31st March 2011

3
Year
Warranty

Meditech ABPM-04 Premium Ambulatory Blood Pressure Monitor

ABPM-04



Price
£1295.00

Was
~~£1389.00~~

- Clinically validated accuracy to BHS protocols and featured on the British Hypertension Society website
- One of the smallest and lightest ABPMs available in the UK
- Proven reliability and ease of use
- Uses popular CardioVisions software, making it easy to programme in minutes and attach reports to popular clinical management systems including EMIS, SystemOne and Vision
- Supplied with adult and large adult cuffs, 2 sets of high capacity rechargeable batteries, battery charger, carry pouch, cable, carry case, a 3 year warranty on monitor and lifetime FREE Premier product support
- Demo video and more information available at www.pmsinstruments.co.uk

InstantCheck Handheld ECG Monitor

InstantCheck



Price
£291.00

Was
~~£330.00~~

- Patient activated ECG recorder for long term monitoring
- Ideal for opportunistic AF screening or loaning to patients
- ECG is recorded by placing thumbs on built in contact electrodes
- 100 thirty second strips can be recorded in memory with date and time
- Simple PC transfer of ECG to PC via USB cable supplied
- JPEG reports can be attached to popular clinical management systems
- Demo video and more information available at www.pmsinstruments.co.uk

Prices exclude carriage and VAT. Offers available until 31st March 2011

ORDER HOTLINE 01628 773233

www.pmsinstruments.co.uk

CardioClip/CardioBlue ECG Event Recorder with Analysis

CardioClip/CardioBlue



Price
£895.00

Was
~~£995.00~~

- Miniature looping pre and post event ECG recorder. Bluetooth version also available
- Stores up to 60 minutes of ECG in digital memory
- Patient activated or programmable event storage
- Ideal for screening patients with intermittent arrhythmias not presenting in clinic
- Uses popular CardioVisions software making it easy to programme and attach reports to popular clinical management systems
- Report software automatically analyses ECG and classifies into normal or abnormal beats and generates a textual report
- Free lifetime Premier product support
- Demo video and more information available at www.pmsinstruments.co.uk

A&D TM-2655P Waiting Room Blood Pressure Monitor

A&D TM-2655



- Clinically validated AA grade validation to BHS protocols
- Patient self measurement solution for routine or opportunistic screening
- Prints results on small paper strip
- Irregular heartbeat feature
- Frees up valuable consultation time
- Energy saving function when not in use
- Supplied with stand, paper, inner cuff covers and security locking cable
- Shown with optional 6121A professional scales at £395 + VAT

Price
£1795.00

Includes 10 rolls of paper
worth £76.00

A&D UA-767+30

Digital Blood Pressure Monitor

UA-767+30



Price
£36.04

Free DT-502EC Digital
Thermometer worth £8.99

- Clinically validated BHS listed BP monitor
- 30 reading memory with average
- Irregular heartbeat indicator
- WHO blood pressure classification indicator
- Supplied with adult 22-32cm cuff, batteries and carry case
- Optional mains adaptor and large cuff 32-45cm available
- 5 year warranty

CardioTens Combined ABPM/ECG Recorder

CardioTens



Price
£1995.00

Was
~~£2295.00~~

- Records 24 hour ABPM profile with sequential recording of ECG strips
- May be used as an ABPM monitor, ECG event recorder or both
- The unit performs continuous real-time beat-to-beat ECG analysis and features automatic event detection with adjustable preset limits and automatic ECG reporting
- Supplied with adult and large adult cuff, CardioVisions software, PC cable, ECG lead, electrodes and 2 year warranty
- Weight approx 240g

PMS
P.M.S. (Instruments) Ltd

Waldeck House, Waldeck Road, Maidenhead, Berkshire SL6 8BR
Telephone: 01628 773233 Facsimile: 01628 770562

Email: sales@pmsinstruments.co.uk
Website: www.pmsinstruments.co.uk



Prices exclude carriage and VAT. Offers available until 31st March 2011

FAX HOTLINE 01628 770562

Email: sales@pmsinstruments.co.uk