

NEW Medset Flashlight PC Based ECG System**Medset Flashlight**

Price	USB version	Bluetooth version
	£1795.00	£2895.00

- Network ready
- Powered by USB port. No batteries. Wireless Bluetooth version also available
- Uses validated German Hanover HES ECG interpretation protocol developed by Cardiologists to provide accurate ECG measurements and interpretation
- Generates high quality PDF or JPEG reports which can be attached to EMIS, SystemOne, INPS clinical management systems for a paperless solution
- Anti drift baseline, muscle and mains power filters for optimum high quality ECGs
- Comparison ECG feature
- Software checks ECG quality before recording
- Record up to 180 seconds of ECG
- 12 channel ECG system with user configurable reports

NEW SB-200 Professional Spot Check Finger Oximeter**SB-200**

- Bright OLED two colour display is easy to read
- Displays SPO2, heart rate and pulse strength
- Changeable display viewing direction
- 2 AAA batteries - last for up to 16 hours continuous use
- Audible alarm function if SPO2 below 90%
- Weighs just 37g size 63.5 x 34 X 35mm

Price
£89.00

Prices exclude carriage and VAT

Meditech ABPM-05

Premium Ambulatory Blood Pressure Monitor

ABPM-05



Price

£1589.00

- Clinically validated accuracy to BHS protocols and featured on the British Hypertension Society website
- One of the smallest and lightest ABPMs available in the UK
- Proven reliability and ease of use
- Uses popular CardioVisions software, making it easy to programme in minutes and attach reports to popular clinical management systems including EMIS, SystemOne and Vision
- Supplied with adult and large adult cuffs, 2 sets of high capacity rechargeable batteries, battery charger, carry pouch, USB cable, carry case, a 2 year warranty and lifetime FREE Premier product support
- Wireless Bluetooth version also available
- Weight 190g size 98 x 69 x 29mm

InstantCheck Handheld ECG Monitor

InstantCheck



Special Offer

£330.00

Save

£50.00

- Patient activated ECG recorder for long term monitoring
- Ideal for opportunistic AF screening or loaning to patients
- ECG is recorded by placing thumbs on built in contact electrodes
- 100 thirty second strips can be recorded in memory with date and time
- Simple PC transfer of ECG to PC via USB cable supplied
- JPEG reports can be attached to popular clinical management systems

Prices exclude carriage and VAT

ORDER HOTLINE 01628 773233

www.pmsinstruments.co.uk

CardioClip/CardioBlue

ECG Event Recorder with Analysis



Special Offer
£995.00

Save
£730.00

- Miniature looping pre and post event ECG recorder. Bluetooth version also available
- Stores up to 60 minutes of ECG in digital memory
- Patient activated or programmable event storage
- Ideal for screening patients with intermittent arrhythmias not presenting in clinic
- Uses popular CardioVisions software making it easy to programme and attach reports to popular clinical management systems
- Report software automatically analyses ECG and classifies into normal or abnormal beats and generates a textual report
- Free lifetime Premier product support

A&D TM-2655P

Waiting Room Blood Pressure Monitor



A&D TM-2655

- Clinically validated AA grade validation to BHS protocols
- Patient self measurement solution for routine or opportunistic screening
- Prints results on small paper strip
- Irregular heartbeat feature
- Frees up valuable consultation time
- Energy saving function when not in use
- Supplied with stand, paper, inner cuff covers and security locking cable
- Shown with optional 6121A professional scales at £395 + VAT

Price
£1795.00

Hygienic Cuff Barriers

Cuff Barriers



- Single patient use cuff barriers for use with any blood pressure cuff providing hygienic protection
- Helps reduce cuff clean up
- Especially suitable for 24 hour ambulatory blood pressure monitors
- Soft tissue against skin ensures comfort
- Latex free
- Call 01628 773233 for free samples
- Paediatric 33 x 13cm sold in a box of 300
- Adult 46 x 17cm sold in a box of 300
- Large adult 57 x 20cm sold in a box of 150

Price	Paediatric per box	Adult per box	Large adult per box
	£49.00	£49.00	£39.00

CardioTens Combined ABPM/ECG Recorder

CardioTens



- Records 24 hour ABPM profile with sequential recording of ECG strips
- May be used as an ABPM monitor, ECG event recorder or both
- The unit performs continuous real-time beat-to-beat ECG analysis and features automatic event detection with adjustable preset limits
- Supplied with adult and large adult cuff, CardioVisions software, PC cable, ECG lead, electrodes and 2 year warranty
- Weight approx 240g

Price	Save
£1995.00	£500.00



Waldeck House, Waldeck Road, Maidenhead, Berkshire SL6 8BR
Telephone: 01628 773233 Facsimile: 01628 770562

Email: sales@pmsinstruments.co.uk
Website: www.pmsinstruments.co.uk



Prices exclude carriage and VAT

FAX HOTLINE 01628 770562

Email: sales@pmsinstruments.co.uk