

Type	Stethoscope with double chest piece, switchable, from brass, chromium-plated and dual membrane, for adults UMDNS Code 13-755 Medical Product Class 4		
Content of delivery	Stethoscope complete with 2 spare ear tips in decorative cardboardbox		
Dimensions of Stethoscope	length (total)		690 (mm)
	length of molded Y-tube		470 (mm)
	diameter of binaurals		5 (mm)
	length of binaurals		190 (mm)
	max. opening of binaurals		300 (mm)
Dimensions of Chest piece	diameter of chest piece bell shaped side		52 (mm)
	height of chest piece		30 (mm)
	diameter of diaphragm		24 (mm)
			45 (mm)
Spring power of binaurals	0,6 mm spring at 150 mm angle width		260 (g)
	at 200 mm angle width		300 (g)
	0,7 mm spring at 150 mm angle width		450 (g)
	at 200 mm angle width		600 (g)
Material	chest piece / binaurals	brass, polished and chromium-plated	
	ear tips super soft	silicone rubber, soft	
	diaphragm, outside	PET, 0.25 mm thick	
	diaphragm inside	CuBe 2, 0.6 mm thick, hard	
	non-chill-rings	NBR	
	molded Y-tube	soft PVC, colour black, latex free, resistant against common detergents and disinfection	
Weight	chest piece		100 (g)
	complete without packing		180 (g)
Dimensions of box	(L x W x H)	370 x 188 x 36	(mm)
Freight data	gross weight (incl. packing) / piece		425 (g)
	volume / piece		2500 (cm ³)
Individual marking	of chest piece maximum 26 letters		
Warranty	5 years		
Manufacturer	ERKA. Kallmeyer Medizintechnik GmbH & Co KG Im Farchet 15 • D 83646 Bad Toelz Tel. +49 (0)8041 8009-0 • Fax: +49 (0)8041 8009-39 • E-Mail: info@erka.org		
Country of origin	Germany		