

PULSE OXIMETER

Health at your Fingertip! SpO2 reads oxygen concentration noninvasively

The Handheld & Fingertip Pulse Oximeter have everything you expect to see in a device used for SpO2 measures a patient's SpO2 and Heart rate.

The unit uses a large display making it easy for the user to read oxygen saturation and Heart rate.

The additional bar graph gives important information about a patient's pulse strength at the probe site.

The unit has a saturation dependent, suppressible pulse tone. Fingertip SpO2 is used for health management & easy application whenever & wherever you want!!



1 **PC-link**
SpO2 and pulse rate PC-printout (Model: SA210 only)

3 **Memory Display**
Could be programmed for 4 intervals including 5/10/30/60-minute "0" stands for no memory setup

5 **LCD Backlight Button**
Adjustable Backlight

7 **Power On/Off Button**
The beep sound will last for 2 seconds to signal the Power On of the device

2 **High / Low SpO2 & Heart Rate**
SET Press ▲ or ▼ button for adjustment

4 **Pulse Strength**
Pulse indicator shows the strength of heartbeat bottom up including the frequency of heartbeat

6 **Sound Off**
The setup of sound volume will be signaled with a beep for 1 second except for the mode of mute

8 **Alert Button**
The alarm will go off while the value exceeds the alert threshold


SA 210
● While SpO2 is low, flashing yellow backlight is on
● While probe is malfunctioning, "Er" signal is displayed



Palpus
Handheld
SA 210

Fingertip
SB100


Fingertip
SB 220



- 1 **Pulse Strength Indicator**
shows the strength of heartbeat
- 2 **%SpO₂ & Heart Rate**
Big LED screen display for easy reading
- 3 **One-touch ON/OFF key**
for easy operation
- 4 **String-attached**
Economical & Portable health management tool

SB100

- ⊙ Reading flashing while SpO₂ is low
- ⊙ While LED or sensor is malfunctioning, displaying "Er" signal



- 1 **Pulse Strength Indicator**
shows the strength of heartbeat
- 2 **%SpO₂ & HR**
2-Way & 2-color OLED display
- 3 **One-touch ON/OFF key**
for easy operation
- 4 **String-attached**
Economical & Portable health management tool
- 5 **Audible alarm system**
Alarm goes off once SpO₂ is 90% below

SB 220

- ⊙ While SpO₂ is low, flashing "exclamation mark" and readings is displayed
- ⊙ While LED or sensor is malfunctioning, "Er" signal is displayed

SPO2 PROBE FAMILY



Adult probe PA100



Child probe PC100



Tube probe PF100



Disposable probe PD100

Rossmax Innotek Corp.

12F, 189, Kang Chien Rd., Taipei, 114, Taiwan.

Tel : +886-2-2659-7888 / Fax : +886-2-2659-7666

E-mail: ling_leonard@rosmaxhealth.com <http://www.rosmaxhealth.com>

CP058220000000021



PULSE OXIMETER

Health at your Fingertip! SpO₂ reads oxygen concentration noninvasively

- Ideal for oxygen therapy assessment
- Ideal for sleep apnea screening
- Ideal for health management
- Ideal for spot-check & continuous monitoring

