

Type	Stethoscope with chest piece from brass, chromium-plated and dual membrane, for adults		
	UMDNS Code	13-755	
	Medical Product	Class 4	
Content of delivery	Stethoscope complete with 2 spare ear tips in decorative cardboardbox		
Dimensions of Stethoscope	length (total)	700 (mm)	
	length of molded Y-tube	470 (mm)	
	diameter of binaurals	6 (mm)	
	length of binaurals	190 (mm)	
	max. opening of binaurals	300 (mm)	
Dimensions of Chest piece	diameter of chest piece	52 (mm)	
	height of chest piece	24 (mm)	
	diameter of diaphragm	45 (mm)	
Spring power of binaurals	0,6 mm spring at 150 mm angle width	260 (g)	
	at 200 mm angle width	300 (g)	
	0,7 mm spring at 150 mm angle width	450 (g)	
	at 200 mm angle width	600 (g)	
Material	chest piece / binaurals	brass, polished and chromium-plated	
	ear tips super soft	silicone rubber, soft	
	diaphragm, outside	PET, 0.25 mm thick	
	diaphragm inside	CuBe 2, 0.6 mm thick, hard	
	non-chill-ring	NBR	
	molded Y-tube	soft PVC, colour black, latex free, resistant against common detergents and disinfection	
Weight	chest piece	105 (g)	
	complete without packing	215 (g)	
Dimensions of box	(L x W x H)	370 x 188 x 36 (mm)	
Freight data	gross weight (incl. packing) / piece	459 (g)	
	volume / piece	2500 (cm ³)	
Individual marking	of chest piece, maximum 20 letters		
Warranty	5 years		
Manufacturer	ERKA. Kallmeyer Medizintechnik GmbH & Co KG Im Farchet 15 • D 83646 Bad Toelz Tel. +49 (0)8041 8009-0 • Fax: +49 (0)8041 8009-39 • E-Mail: info@erka.org		
Country of origin	Germany		