




PMS (Instruments) ltd
 Waldeck House, Waldeck Rd, SL6 8BR
 Tel. 01628 773233
 Fax. 01628 770562

Patient data			
name	Normal CardioTens		
ID	4298-2349	birth date	22/11/1959
Mother name			
Birth place		Sex	Male





Blood pressure list						
Time	T	SBP mmHg	DBP mmHg	pulse bpm	remark	note
09:35	M	136	095	070		
09:40	X	117	060	118		
10:00	A	133	075	073		
10:20	A	125	083	060		
10:40	A	123	082	064		
11:00	A	125	088	071		
11:20	A	128	085	074		
11:40	A	124	088	073		
12:00	A	128	089	081		
12:20	A	131	086	072		
12:40	A	125	071	079		
13:00	X	149	123	200		
13:20	A	118	081	065		
13:40	A	116	085	065		
14:00	A	120	082	074		
14:20	A	117	077	082		
14:40	A	122	087	088		
15:00	A	119	085	084		
15:20	-					no data
15:40	-					no data
16:00	-					no data
16:20	?					cuff leaking or loose
16:40	A	128	081	078		
17:00	A	132	077	078		
17:20	A	110	085	081		


Blood pressure list						
Time	T	SBP mmHg	DBP mmHg	pulse bpm	remark	note
17:40	A	130	082	070		
18:00	A	120	085	097		
18:20	A	113	086	082		
18:40	A	126	091	088		
19:00	A	125	082	089		
19:20	A	118	070	073		
19:40	A	115	074	079		
20:00	A	117	075	074		
20:20	A	117	076	083		
20:40	A	108	066	064		
21:00	A	120	084	083		
21:20	A	117	088	087		
21:40	A	122	090	079		
22:00	A	149	111	096		
22:30	A	116	085	071		
23:00	A	124	084	075		
23:30	A	120	085	060		
00:00	A	122	083	061		
00:30	A	108	073	049		
00:39	E					event button pressed
01:00	A	101	072	049		
01:30	A	111	076	051		
02:00	A	102	071	051		
02:30	A	101	066	060		
03:00	A	112	071	061		
03:30	A	109	064	049		
04:00	A	097	063	050		
04:30	A	104	057	053		
05:00	A	101	053	049		
05:30	?					cuff leaking or loose
06:00	?					failed reading
06:15	A	102	063	049		
06:30	A	121	075	065		
06:45	A	096	060	060		
07:00	A	089	051	055		
07:15	A	100	064	056		
07:30	A	110	069	063		
07:45	A	106	073	080		
08:00	A	129	076	107		
08:20	A	119	079	082		
08:40	A	122	069	082		
09:00	A	128	089	075		
09:20	A	127	091	092		
09:40	-					no data
14:45	-					no data
15:00	-					no data
15:15	-					no data
15:30	-					no data
15:45	-					no data
16:00	-					no data
16:15	-					no data
16:30	-					no data

Blood pressure list						
Time	T	SBP mmHg	DBP mmHg	pulse bpm	remark	note
16:45	-					no data
17:00	-					no data
17:15	-					no data
17:30	-					no data
17:45	-					no data

ABPM summary (09/02/2005 09:40 - 10/02/2005 09:40)				
period	Overall	Morning	Day	Night
time range	Overall	06:00 - 08:00	08:00 - 22:00	22:00 - 06:00
SBP/DBP high limit	~127/77~	135/85	135/85	120/70
SBP/DBP weighted average	116/77	103/65	122/82	111/74
Pulse weighted average	70	61	78	58
readings	58	7	36	15
Successfully ratio %	88	78	86	100
SBP/DBP D/N index %	9/10			
morning surge	3			



ECG spikechart



Holter summary			
Channels / sampling rate	2 / 200 Hz	Total beats	18252
		Of which:	
The registration started	09/02/2005 09:35	Normal	17807 (97.6%)
The registration ended	10/02/2005 03:04	VES	87 (0.5%)
Total registration time	03 h 37 min	SVES	358 (2.0%)
Valuable length	03 h 37 min	PM	0 (0.0%)
Total lead off / noise	00 h 00 min	Other	0 (0.0%)
Pacemaker spike detection	disabled	Artifact	0 (0.0%)

Heart rates, RR intervals	Average	Max	Min	Max time	Min time
HR [1/min]	77	125	45	09/02/2005 15:18	10/02/2005 00:59
HR equivalent intervals [ms]	779				
RR filtered [ms]	727	1730	415	10/02/2005 02:33	09/02/2005 15:21
RR unfiltered [ms]	727	1730	415	10/02/2005 02:32	09/02/2005 15:21

Tachy / brady / pause	Number of events	Longest event	
		Time	Length
Tachycardia	164	09/02/2005 15:56	2 min 57 s
Bradycardia	65	10/02/2005 00:56	1 min 56 s
Pause (2.5 s)	0		

QRS statistics						
	Total	Max/h	Longest		Highest rate	
			Duration	HR [1/min]	Duration	HR [1/min]
Ves						
Isolated	87	39				
Couplet	0	0				
Triplet	0	0				
V-tachycardia	0	0				
Bigeminy	12	8	3 s	94	1 s	121
Trigeminy	7	3	5 s	63	2 s	121

QRS statistics						
SVES						
Isolated	12	3				
Couplet	1	1				
Triplet	2	1				
SV-tachycardia	9	3	00:00:24	114	00:00:24	114

