

ERKA. Perfect-Aneroid

Made in Germany

CE 0124



with screw valve for sensitive adjustment

scale 48 mm Ø, in case

with velcro cuff
art. no. 201 003 02

with velcro cuff washable
art. no. 201 004 02

with hook cuff
art. no. 201 001 02

scale 56 mm Ø, in case

with velcro cuff
art. no. 203 003 02

with hook cuff
art. no. 203 001 02

Our features

- Manometer housing brass chrome-plated
- Membrane-protecting 2-tube technology
- Measuring accuracy over the total range: $\pm 3\text{mmHg}$

with patented automatic valve

scale 48 mm Ø, in case

with velcro cuff
art. no. 202 003 02

with velcro cuff washable
art. no. 202 004 02

with hook cuff
art. no. 202 001 02

scale 56 mm Ø, in case

with velcro cuff
art. no. 204 003 02

with hook cuff
art. no. 204 001 02

Your convenience

- Flat, easy-to-handle precision manometer, ideal for the doctor's bag
- All materials recyclable

Clinic model with protective cover

scale 48 mm Ø, without case

with velcro cuff washable
art. no. 215 004 02

Further variants

Perfect Aneroid with screw valve
scale 48 mm Ø,
in case with velcro cuff

for children
art. no. 201 063 02

for infants
art. no. 201 073 01

for babies
art. no. 201 083 01

for thighs
art. no. 201 043 02

Subject to technical alterations

Rev. Std. 004 / Jan. 1998