



**ERKA. Profi**  
Made in Germany  
**CE** 0124

#### Our features

- Precision manometer, Aluminium, black anodised  
Scale-Ø 48 or 56 mm with sensitive air-reduction valve
- Scale, 300 mmHg  
2 mm scale divisions
- New technology with air intake via the manometer capsule in the specially tempered copper-beryllium membrane
- Black, kind to the skin velcro cuff made of 100 % cotton suitable for washing or dry cleaning

#### Your convenience

- Classical model for the professional
- Elegant black anodised finish
- Sensitively adjusted air-pressure relief valve, high-quality chrome finish
- User-friendly. The ergonomic design facilitates visually clear easy reading of the blood pressure values
- Convincing ERKA quality

#### For your Order

Profi  
Sphygmomanometer  
black anodised finish  
with black velcro cuff  
twin-tube model

available in

Ø 48 mm Art. No. 228 053 02  
Ø 56 mm Art. No. 229 053 02

Subject to technical alterations

Rev. Std. 002 / January 1998

**ERKA.**

***P.M.S (Instruments) Ltd***

Waldeck House, Waldeck Road  
Maidenhead, Berkshire, SL6 8BR  
Tel: 01628 773233 or Fax: 770562  
E-Mail: Sales@pmsinstruments.co.uk