

UA-1000 SERIES

UA-1010 / UA-1020 / UA-1030T



Blood Pressure Monitors



Advanced Cuff design
Advanced Pressurization
Advanced Deflation



JAPAN
Quality & Design



A&D

A&D Medical
<http://www.aandd.jp>

Technology You Can Trust



UA-1000 SERIES

UA-1010 / UA-1020 / UA-1030T

A&D, a pioneer in home blood pressure monitoring, has once again come up with a new series of world-class monitors with 3A Technology that sets new industry standards and enhances value for users.

SlimFit cuff



UA-1010



Simple monitoring with more accuracy It offers easy, one-button operation with new functions that will guide user to accurate measurements. UA-1010 uses the new latex-free and metal-free SlimFit cuff which is more friendly to sensitive skin.

UA-1020



Individualized monitoring with more comfort UA-1020 is equipped with newly designed SmoothFit cuff, which enables painless, comfortable measurement. Also it has two unique modes. MyPressure mode provides individualized pressurization and shorten the measurement time. The TriCheck™ mode automatically takes three consecutive readings with one-minute intervals for most accurate results.



SmoothFit cuff
*Patent pending



SmoothFit cuff
*Patent pending

UA-1030T



Talking monitor for everyone UA-1030T has come with the latest talking technology with all features of UA-1020. It announces the results in 5 languages (English, French, Spanish, Italian, and Chinese). You can increase volume through six levels or switch on to silent mode if required.



Cuff Fit Error Indicator



Movement Error Indicator



Irregular Heartbeat Indicator (IHB)



Latex-free and metal-free, SlimFit Cuff



WHO Blood Pressure Classification Indicator



M 90 90 memories + Average Reading (All/AM/PM)



Clock & Calendar



Newly designed SmoothFit Cuff for painless and comfortable measurement (Latex-free and metal-free)



My Pressure Setting mode (AUTO/180/210/240)



TriCheck™ mode for consecutive 3 measurements and average reading



Clear voice announcement of measurement results



Advanced Cuff design
Advanced Pressurization
Advanced Deflation

Latex-free and metal-free plus wrinkle-free cuff designs.

User selection made possible between A&D's fuzzy logic inflation and anxiety-free "MyPressure" setting.

Linear-motor driven active deflation valve makes it measure accurately for all cuff sizes.



A&D Company, Limited

3-23-14 Higashi-Ikebukuro, Toshima-ku, Tokyo 170-0013 Japan
Telephone: [81](3)5391-6132 Fax: [81](3)5391-6148
1-243, Asahi, Kitamoto-shi, Saitama, 364-8585 Japan
Telephone: [81](0485)93-1111 Fax: [81](0485)93-1119
(EN46001 Certified Factory)
<http://www.aandd.jp>

A&D INSTRUMENTS LTD.

P.M.S (Instruments) Ltd

Waldeck House

Waldeck Road

Maidenhead

SL6 8BR

Tel 01628 773233 Fax 01628 770562

Email sales@pmsinstruments.co.uk