

DIGITAL BLOOD PRESSURE MONITOR

UA-767

Automatic

DIGITAL BLOOD PRESSURE MONITOR UA-767

Automatic

One Button Operation

The UA-767 has one button operation : Power-on, start, pressurization at pre-set level, automatic measurement, display of systolic and diastolic blood pressure readings and heart rate, automatic cuff deflation and unit shut off are all accomplished with the push of only one button.

Easy to Read, Liquid Crystal Display


Like its brother model, UA-777, this unique display simultaneously displays systolic /diastolic blood pressure and heart rate, giving the user instantaneous readings of all critical data

Advanced Oscillometric Measurement Technique

The same Oscillometric technique used by physicians is incorporated into the UA-767. This technique, along with advanced software and an extra sensitive pressure sensor, assures reliable and accurate blood pressure measurement.

Ultra Modern Design

The UA-767 is the most attractive digital blood pressure monitor on the market. It can be used with pride in either home or office.



Advanced ONE BUTTON OPERATION

- One button operation : automatic power on cuff inflation blood pressure measurement and cuff deflation.
- Easy to read, liquid crystal display : Systolic, diastolic and heart rate shown simultaneously.
- Advanced Oscillometric method for hospital accuracy.



AND

A&D MEDICAL

Clearly a Better Value

**Call P.M.S (Instruments) on
01628 773233 for sales.**